

COMBINED DEVICE MALFUNCTION & SAFETY INCIDENT REPORT

PATIENT: Female | Age: 62 | ID: SYN-9005

PRODUCT: Rheumaject 50mg Auto-Injector | Lot: RJ-99381

PROBLEM: Needle shield stuck during self-administration.

REACTION: Severe subcutaneous hematoma at injection site.

SEVERITY: Serious = Yes | Hospitalization = No | Life-Threatening = No

NARRATIVE: Mechanical device failure caused needle bending and severe tissue trauma.